



Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 1 of 31	Standard Operating Procedure	

***INVITTOX* Protocol**

epiCS[®]

**Skin Sensitisation and Potency Test
Method**

(SSPT method)

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 2 of 31	Standard Operating Procedure	

A. Protocol Introduction

The format of this INVITTOX protocol was adapted from the ECVAM Database Service on Alternative Methods to Animal Experimentation. The content of the protocol is based on a SOP kindly provided by Prof. Dr. Sue Gibbs (The Netherlands) and Prof. Dr. Emanuela Corsini (Italy) "Identification and ranking potency of a contact sensitiser using the *in vitro* skin equivalent epiCS". The protocol is established and approved by CellSystems GmbH.

epiCS[®] Skin Sensitisation and PotencyTest (SSPT)

The epiCS SSPT method is used to identify contact allergens and to rank them according to their potency (extreme, strong, moderate and weak sensitising potency) with the aid of the human skin equivalent epiCS (CellSystems) and the release of interleukin-18 (IL-18) into the culture medium. The test system is available as a kit, comprising reconstructed epidermis, culture media and culture plates. The IL-18 ELISA kit is not included.

Please, contact the CellSystems experts before you start this test, to help you establish it at your site.

OBJECTIVES & APPLICATIONS

Type of Testing:

Replacement, as standalone test or within an Integrated Testing Strategy

Level of Toxicity Assessment:


Hazard identification, toxic potential, contact allergen identification and potency determination

Purpose of Testing:

Classification and labelling, safety. Determine the allergenic potential and potency (strength) of low molecular weight substances.

Context of Use:

Currently no validated *in vitro* method exists to replace animal testing for the identification of skin sensitising substances. Several cell-based methods are under development and might be required to be used in a test battery to fully replace animal testing. To date, the mouse LLNA is the preferred and accepted method for assessing skin sensitisation of most substances.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 3 of 31	Standard Operating Procedure	

Applicability Domain:

The test described in this protocol is designed for the classification of substances and was established for liquids, viscous and solid test substances.

BASIS OF THE METHOD

Skin sensitisation is an immune mediated reaction to low molecular weight substances, involving several cell types which actively participate in contact allergy within the human skin (Cavani 2007). Keratinocytes, the major cell type in human epidermis, play a role in all phases of allergic contact dermatitis, from the early initiation phase with penetration through the stratum corneum and the secretion of inflammatory cytokines which are required for Langerhans' cells migration and T-cell trafficking, leading to the development of allergic contact dermatitis. In vivo, the local lymph node assay (LLNA) is used to study the induction phase of skin sensitisation and provides quantitative data suitable for dose-response assessment (OECD TG 429, 2010).

This protocol of in vitro skin sensitisation test using epiCS may also enable to identify contact allergens and to rank them according to their sensitising potency (extreme, strong, moderate and weak). The advantage of the assay is the possibility to identify in the same test the allergenic potential of a substance and its potency. While topically applied, it mimics human exposure and overcomes all drawbacks of traditional submerged culture, including substance solubility and stability in culture medium.

The test consists of a topical exposure of test substances to a human reconstituted epidermis model followed by a cell viability test to determine the EC₅₀ value. Substances can be diluted in the solvent that shows best solubility. Neat mixtures and finished products might be tested for their allergenic potential without dilution prior to experimentation.


Cell viability is measured by dehydrogenase conversion of MTT [(3-4,5-dimethyl thiazole 2-yl) 2,5-diphenyl tetrazolium bromide], present in cell mitochondria, into a blue formazan salt that is quantitatively measured after extraction from tissues. The reduction of viability of tissues exposed to substances in comparison to negative controls (treated with solvent) is used to determine to calculate EC₅₀.

Interleukin-18 (IL-18) is constitutively secreted by keratinocytes and may be induced by test substances to different degrees. The IL-18 production level induced by a test substance is one parameter of the prediction model for this assay.

EXPERIMENTAL DESCRIPTION

Endpoint & Endpoint Detection:

Cytotoxicity is measured by MTT assay and the release of IL-18 into the culture medium by ELISA.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 4 of 31	Standard Operating Procedure	

MTT Assay: Cell viability is used as endpoint and measured by dehydrogenase conversion of MTT [(3-4,5-dimethyl thiazole 2-yl) 2,5-diphenyltetrazoliumbromide], present in cell mitochondria, into a blue formazan salt that is quantitatively measured after extraction from tissues. The reduction of viability of tissues exposed to substances in comparison to negative controls (treated with solvent) is used to calculate EC₅₀.

IL-18 ELISA: A commercially available IL-18 ELISA kit (#7620, MBL International Corporation, Japan) can be used to determine the secretion of the cytokine into the culture medium within 24 hours.

Endpoint Value:

Results from the dose response experiments are expressed as fold increase in IL-18 release in order to determine the allergenic potential. A substance is classified as contact allergen if it induces ≥ 5 fold increase in IL-18 release (ELISA) compared to solvent treated control at cell viability $\geq 5\%$ and $\leq 40\%$ (EC₅₋₄₀) in at least one of the concentrations tested.


For sensitiser potency, the EC₅₀ value and IL-18 (SI-2) value is determined. The EC₅₀ value and IL-18 (SI2) of the unknown substance is correlated to the values obtained for a standard test panel of substances (dos Santos et al., 2011).

Test System:

epiCS (CellSystems GmbH, Troisdorf, Germany) consists of normal, human-derived epidermal keratinocytes, which have been cultured to form a multilayered, highly differentiated model of the human epidermis. It consists of organized basal, spinous and granular layers, and a multilayered stratum corneum containing intercellular lamellar lipid layers *arranged* in patterns analogous to those found in vivo. The epiCS tissues (surface 0.6 cm²) are cultured on specially prepared cell culture inserts and shipped world-wide as kits, containing tissues on shipping agarose together with epiCS Culture Medium, and 6-well plates.

Basic Procedure:

On day of receipt, epiCS are conditioned by incubation for release of transport-stress related compounds and *debris* at least 4 h or overnight. After pre-incubation, tissues are supplied with new pre-warmed culture medium. Tissues are topically exposed to the test substances impregnated in 8 mm diameter paper filters (designed for the use with 8 mm Finn Chambers). The paper filters were loaded with 100 μ l for water soluble compounds or 25 μ l for ethanol or AOO soluble compounds) on a glass plate under sterile conditions, excess liquid removed, and applied on to the tissues surface for 24 hours. Hereafter, cytotoxicity is measured by MTT assay and the release of IL-18 into culture medium by ELISA. Results from the dose response experiments are expressed as fold increase in IL-18 release in order to determine the allergenic potential.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 5 of 31	Standard Operating Procedure	

MTT assay: Transfer the tissues to 24 well plates containing 300 µl MTT medium (MTT 1 mg/ml). After a 3 hr MTT incubation the blue formazan salt formed by cellular mitochondria is extracted for 2h (RT) or overnight sealed (2-8°C) with 2 ml isopropanol per tissue and the optical density of the extracted formazan is determined in a spectrophotometer at 540-570 nm. Relative cell viability is calculated for each tissue as % of the mean of the negative control tissues.

Correction for MTT reducers: Please follow the test for interference of substances with MTT as described in the epiCS Skin Irritation INVITTOX Protocol (see Annex 1).

IL-18 ELISA: Please follow the instructions of the IL-18 ELISA kit (#7620, MBL International Corporation, Japan).

DATA ANALYSIS / PREDICTION MODEL

The test protocol allows predicting the skin sensitisation potential of unknown test substances according to the following prediction model:

Allergenic potential:

≥ 5 fold increase in IL-18 release (ELISA) compared to solvent treated epiCS at cell viability ≥ 5 % and ≤ 40 % (EC₅₋₄₀).


Potency assessment:

Primary parameter: Cytotoxicity (MTT assay) expressed as EC₅₀ value (EC₅₀ = effective substance concentration required to reduce epiCS metabolic activity - corresponding to cell viability - to 50% of the maximum value compared to solvent exposed tissues).

Secondary parameter: IL-18 Substance concentration resulting in : ≥ 2 fold increase in IL-18 release (ELISA) compared to solvent treated tissues (SI-2).

TEST COMPOUNDS & RESULTS SUMMARY

The assay is performed according to the SOP in order to determine the EC₅₀ value of a substance of (un)known identity, solubility and molecular mass, and the release of IL-18 associated with its treatment. First the maximum solubility of the substance is identified by dissolving the compound in both AOO and in 1% DMSO in epiCS Culture Medium. Dilutions are made in the following order: 200 mg/ml, 150 mg/ml, 100 mg/ml, 50 mg/ml, 25 mg/ml, 12.5 mg/ml, 6.25 mg/ml, 3.125 mg/ml, 1.56 mg/ml, 0.78 mg/ml, 0.39 mg/ml, 0.20 mg/ml, 0.10 mg/ml until a clear solution is reached. The solvent with the highest dissolving capacity (AOO (4+1) or 1% DMSO) is chosen. Substances are tested in a dose response starting with the highest soluble concentration (preferably max. 200 mg/ml) and decreasing with 2-fold serial dilutions until 0.1 mg/ml (or the closest lower concentration) is reached. Whereas substance dose responses are tested in single fold in each independent experiment, control conditions (unexposed, solvent(s) and positive control), should preferably be tested in duplicate per independent experiment.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 6 of 31	Standard Operating Procedure	

During the development of this assay 11 coded sensitisers and 9 coded non-sensitisers (11 *in vivo* contact allergens and 9 *in vivo* irritant substances) were tested and the results published in Gibbs et al (2013). Five sensitisers* and 5 non-sensitisers** of this panel were taken from the reference substances list mentioned in the OECD TG 429 “Skin Sensitisation: Local Lymph Node Assay”. A ring trial in three different laboratories is performed with a small subset of substances to evaluate the prediction model (see above) followed by an independent statistical analysis in 2014.

MODIFICATIONS OF THE METHOD

None

DISCUSSION

The epiCS Sensitisation Test is an easy to perform *in vitro* test which only needs minimal training. Besides a standard cell culture and substances laboratory no special equipment is required.


The method offers a high reproducibility due to standardised materials and processes during the production. Furthermore the epiCS is based on human cells to predict effects on humans and may serve as a complete replacement of the *in vivo* LLNA test in mice. In view of stimulated discussions on Integrated Testing Strategies to replace the LLNA assay the epiCS Skin Sensitisation Test might support those testing strategies and thereby complement the assay portfolio to reach this goal.

Limitations: The solvents that have been successfully tested on epiCS are water, PBS, acetone:olive oil (AOO) (4 parts acetone : 1 part olive oil), absolute ethanol and DMSO (maximum 1%). Water insoluble substances can be tested although care should be taken that solvents are not used at irritant concentrations. The substance needs to be in solution in order to be tested.

STATUS

In Development:

The epiCS SSPT method was published by Gibbs et al. (2013).

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 7 of 31	Standard Operating Procedure	

Known Laboratory Use:

epiCS is used in industry and academia for research, efficacy testing and regulatory toxicology testing (skin irritation testing, skin corrosivity testing, sub-categorisation of corrosives)

Participation in Validation Studies:

This protocol was refined by CellSystems GmbH and might be part of a blind trial multicenter validation study. Modifications of the protocol are possible.

Regulatory Acceptance:


To date there is no regulatory accepted *in vitro* skin sensitisation test method available.

PROPRIETARY ISSUES

The intellectual property rights for the test system are held by CellSystems GmbH. No IPRs are associated with the present method.

ABBREVIATIONS & DEFINITIONS

MTT:	3-(4,5-dimethyl-thiazoyl-2-yl)-2,5-diphenyl-tetrazolium bromide
NC:	negative control
OD:	optical density
PC:	positive control
ref.:	reference
RH:	relative humidity
RT:	room temperature
SD:	standard deviation
UN-GHS:	United Nations Globally Harmonized System.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 8 of 31	Standard Operating Procedure	

* Sensitisers from the reference substances list (OECD TG 429):

DNCB
Isoeugenol
2-Mercaptobenzothiazole
Eugenol
Citral

** Non-sensitisers from the reference substances list (OECD TG 429):

Sodium Lauryl Sulfate
Lactic Acid
Methyl Salicylate
Salicylic Acid
Chlorobenzene

B. Technical Description


Procedure Details, Version 1.2, September 2015

epiCS SSPT

Contact Person: **Dr. Oliver Engelking**

Full Address: CellSystems Biotechnologie Vertrieb GmbH
 Langeler Ring 5
 53842 Troisdorf
 Germany

phone: +49 (0) 2241-25515-0
fax: +49 (0) 2241-25515-30
e-mail: oliver.engelking@cellsystems.de

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	
Page 9 of 31	Standard Operating Procedure	

HEALTH & SAFETY ISSUES

All biological components of the epidermis and the culture medium are tested by manufacturer for viral, bacterial, fungal and mycoplasma contamination.

The epidermal cells are taken from healthy volunteers negative for HIV, and Hepatitis B+C virus. Nevertheless, handling procedures for biological materials should be followed:

- a) Wear gloves during handling with the skin and kit components
- b) After use, the epidermis, the material and all media in contact with it should be decontaminated prior to disposal (e.g. using 10 % bleach, special containers or autoclaving).
- c) Examine all kit components for integrity and deviation from standard appearance.

Please contact CellSystems if you have any questions or concerns.

MATERIALS AND PREPARATIONS


TEST SYSTEM

Standard epiCS Kit Components

Sealed 24-well plate (epiCS)	<i>Contains up to 24 inserts with epidermis tissues (origin human keratinocytes) on agarose</i>
6-well plates (sterile)	<i>For assay / culture</i>
epiCS Culture Medium	<i>For tissue culture</i>

Expiration and Kit Storage

reference	description	conditions	shelf life
CS-1001	epiCS (reconstructed epidermis)	refrigerator (2-8°C)	72 h

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	
Page 10 of 31	Standard Operating Procedure	

	epiCS Culture Medium:		
CS-3050	50 ml	refrigerator (2-8°C)	(see label)
CS-3051	75 ml		
CS-3052	100 ml		
CS-3053	125 ml		

Note: Examine all kit components for integrity. If there is a concern call CellSystems immediately.

Phone: +49 (0)2241-25515-0


EQUIPMENT

Fixed Equipment:

Laminar flow hood	<i>For safe work under sterile conditions</i>
Humidified incubator (37±1°C, 5±1 % CO ₂ , 90±10 % relative humidity (RH))	<i>For incubating tissues prior to and during assays</i>
Laboratory balance	<i>For pipette verification and checking substance weight</i>
Vacuum source/trap (optional)	<i>For aspirating media and solutions</i>
96-well plate spectrophotometer	<i>For reading OD</i>
6-well plate and 24-well plate	<i>For culture and MTT assay</i>
Plate shaker	<i>For extraction of formazan</i>
Timers	<i>To be used during application of test materials</i>
Sterile, forceps or tweezers	<i>For handling tissue inserts</i>
37±1°C water bath	<i>For warming Media and MTT solution</i>
Adjustable pipette	<i>For pipetting MTT isopropanol</i>

Consumables:

Sterile paper towel	For blotting of cell culture inserts
Adhesive tape or Parafilm	Covering plates during formazan extraction

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	
Page 11 of 31	Standard Operating Procedure	

96-well plates	<i>For reading OD</i>
6-well plate and 24-well plate	<i>For culture and MTT assay</i>
Finn Chamber filter paper discs 8 mm (Smart Practice; Ref. AL5112GN); website with local distributors: http://finnchambers.com/en/Distributors.aspx	<i>For substance application</i>

MEDIA, REAGENTS, SERA, OTHERS

PBS Ca/Mg-free (sterile)	<i>For MTT-Assay</i>
MTT - Thiazolyl Blue Tetrazolium Bromide (Sigma, ref. M-5655, cell culture tested, purity min. 97,5 %)	<i>For MTT assay</i>
Isopropanol	<i>For extraction of the formazan crystals</i>
DMSO, Acetone, Olive oil, Ethanol	<i>For solvent preparation</i>

PREPARATIONS


Media and Endpoint Assay Solutions:

MTT solution (prepare freshly on day of testing)

Per 24 tissues, dissolve 12.5 mg MTT (Sigma ref. M5655) in 12.5 ml PBS and thoroughly vortex this working solution. Keep the MTT solution in the dark and warm to 37°C prior to usage (use within 2 hrs).

Alternatively, a stock solution might be used as follows:

- Prepare MTT stock solution in PBS (5 mg/ml) and store aliquots at -20°C until use.
- For 50 ml stock solution, weigh 250 mg of MTT and add 50 ml PBS.
- PBS should NOT contain Ca²⁺ (calcium) or Mg²⁺ (magnesium).
- MTT working solution is 1 mg/ml diluted in PBS (2.5 ml stock solution plus 10 ml PBS).

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	
Page 12 of 31	Standard Operating Procedure	

Safety precaution: MTT is toxic (Risk phrases: R26, R68, R22, R36, R37, R38).

Wear protective gloves during manipulation with MTT solution.

ELISA kit

Please follow the instructions of the MBL ELISA kit #7620 (MBL, Japan) with the following exception:

The analysis of the collected culture medium for the IL-18 content should be done **without** dilution (1 insert/well in a six well plate, 1 ml medium for 24h). Only, if OD values are out of the photometer's linearity range, samples should be diluted.

Prepare the standard curve at a concentration range from 0 to 500 pg/ml ideally with equidistant concentrations.

Test Compounds:

Follow safety instructions for respective compounds.


The use of solvent for application onto the tissues differs depending on substance solubility; see chapter "test compounds and results summary".

Positive Control Solution(s):

Dinitrochlorbenzol (DNCB) dissolved in AOO should serve as a positive control. AOO is prepared by combining 4 parts acetone with 1 part of olive oil. The DNCB concentration that will serve as PC needs to be determined in previous experiments. This can be done with the use of the SSPT training kit (CS-1018-SSPT) or when testing in the range of: 0.156 to 10 mg/ml).

Negative Control:

Unexposed epiCS (no paper filter and no solvent) and solvent treated epiCS (with paper filter) serve as a negative control.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 13 of 31	Standard Operating Procedure	

METHOD

TEST SYSTEM PROCUREMENT:

CellSystems Biotechnologie Vertrieb GmbH

Langeler Ring 5

53842 Troisdorf

Germany

phone: +49 (0) 2241-25515-0

fax: +49 (0) 2241-25515-30

e-mail: info@cellsystems.de

ROUTINE CULTURE PROCEDURE:

Day prior to testing:

1. Upon receipt of the shipment, examine all kit components for integrity. If there is a concern call CellSystems immediately.

Contact person:


Dr. Oliver Engelking

CellSystems Biotechnologie Vertrieb GmbH

phone: +49 (0) 2241-25515-0

e-mail: oliver.engelking@cellsystems.de

2. Keep epiCS Culture Medium at 2-8°C.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 14 of 31	Standard Operating Procedure	

Tissue conditioning:

1. *Do not pre-warm epiCS Culture Medium!*
2. Pipette 1 ml of the cool epiCS Culture Medium into the wells of sterile 6- well plates, which will be used during the test (see annex 3)
3. Remove the shipped multi well plate from the package. Under sterile conditions carefully (using sterile tweezers) take out each insert containing the epiCS tissue. Remove any remaining agarose that adheres to the outer sides of the insert by gentle blotting on the sterile blotting paper, and place the tissues in one well of the prepared 6-well plate. Avoid air bubble formation underneath the tissue culture insert.
4. Perform visual inspection of the inserts within next 5 min. Record any tissue defects and excess of moisture on the surface. Do not use defect tissues or tissues with excessive moisture on the surface.
5. Incubate the plates at least for 4 h or overnight 37±1°C, 5±1 % CO₂, 90±10 % RH).
6. Store remaining epiCS Medium at 2-8°C.
7. Prepare and sterilize all devices which will be used in the assay.


TEST MATERIAL EXPOSURE PROCEDURES:

Performing substance exposure experiments involves a number of steps. Each step is addressed individually.

Preparation of solvents

Six different solvents have been tested topically on epiCS.

- AOO (4+1) - 8 ml acetone + 2 ml olive oil
- AOO (1+19) - 0.5 ml acetone + 9.5 ml olive oil
- 1% DMSO - 50 µl DMSO + 4950 µl epiCS culture medium
- PBS – 100% PBS
- Water – sterile dH₂O
- Absolute ethanol

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 15 of 31	Standard Operating Procedure	

Dissolving substances

When dissolving substances in the appropriate solvent, the substances themselves contribute to the total volume of a substance solution by having “their own volume” in the substance solution. This should be taken into account when preparing the substance solution. Therefore, the following applies when handling liquid, solid and viscous substances:

Example liquid substances:

Assumption: 1 µl of liquid substance equals 1 µl in volume in the substance solution. Maximum test concentration of the liquid substance is 200 mg/ml, the solvent is AOO (4+1).

- e.g. you take 75 µl liquid substance → this weighs 65 mg
- to obtain the max. test concentration of 200 mg/ml you need to dissolve the 65 mg in a total volume of 325 µl AOO (calculation $(65 \cdot 1000)/200 = 325 \mu\text{l AOO}$)
- you need to add $325 \mu\text{l} - 75 \mu\text{l} = 250 \mu\text{l AOO}$ to the 75 µl liquid substance.

Example solid and viscous substances:

Assumption: 1 mg solid or viscous substance equals 1 µl in volume in the substance solution.


Max. test concentration of the solid or viscous substance is 200 mg/ml, the solvent is AOO (4+1).

- e.g. you take 92 mg solid or viscous substance → these 92 mg equal 92 µl
- to obtain the max. test concentration of 200 mg/ml you need to dissolve the 92 mg in a total volume of 460 µl AOO (calculation → $(92 \cdot 1000)/200 = 460 \mu\text{l AOO}$)
- you need to add $460 \mu\text{l} - 92 \mu\text{l} = 368 \mu\text{l AOO}$ to the 92 mg solid or viscous substance.

Be aware: When a liquid compound is viscous, handle it as a solid compound. This is more accurate.

Determination of Master Stock Concentration: Maximum solubility of the substance

The maximum solubility of the substance is identified by dissolving the compound in one of the above mentioned solvents. Wherever possible AOO (4+1) should be the preferred solvent of choice. Dilutions are made in the following order: 200 mg/ml, 150 mg/ml, 100 mg/ml, 50 mg/ml, 25 mg/ml, 12.5 mg/ml, 6.25 mg/ml, 3.125 mg/ml, 1.56 mg/ml, 0.78 mg/ml,

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 16 of 31	Standard Operating Procedure	

0.39 mg/ml, 0.20, 0.10 mg/ml until a clear solution is reached. The solvent with the highest dissolving capacity is chosen for the substance application.

This highest test concentration is approved if it does **not** result in

- i) Two layers
- ii) A milky suspension
- iii) Visible precipitation of the substance (after 10 min waiting time at RT)

If any of these effects occur then substance concentrations have to be lowered in the solvent(s) until these effects do not occur anymore. That will be then the highest possible test concentration in the exposure experiments.


General remarks

- Be aware that the volume of substances is taken into account when dissolving the substances
- Substances must be handled with care as they may be harmful
- Be aware that substances may be quite vaporous. Work in appropriate ventilated environment (e.g. acid cabinet or flow cabinet with external air outlet)
- Also make sure that laboratory equipment used during the substance exposure experiment and the analysis after the exposure is properly cleaned and waste is properly removed
 - it is therefore recommended to perform the analysis after the 24 h exposure as well in the laminar flow
- Substances must be freshly made just before application
- The use of gloves, lab-coat and mask is recommended
- For liquid substances, all crystals should be dissolved before pipetting, by leaving the liquid substances at RT for up to 1 hour after removal from cooling (if cooling conditions are required for the storage of the substance).

Preparation of substance test concentrations

N.B. It is known that AOO dissolves some plastic lab devices. It is therefore recommended to use polypropylene tubes to prepare the substance solutions dissolved in AOO and impregnate the paper filters for the substance exposure on devices composed of glass (e.g. glass slides or tissue culture dishes made of glass).

If AOO is the solvent of choice for the substance to be tested, take into account that acetone is a volatile liquid. Therefore, open and close tubes with AOO (+/- a substance) rather quickly to avoid too much evaporation of the AOO. Do not use small volumes < 300 µl. If you do so this will increase error margins.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 17 of 31	Standard Operating Procedure	

- Determination of interference of substance with MTT assay (an alternative method is to follow what is described in the SOP epiCS Skin Irritation Test method , see Annex 1)
- Some substances may interfere with the MTT assay, which is used to determine the EC₅₀ value and therefore needed to be excluded from the epiCS potency assay – these substances fall outside of the applicability domain of the assay. In order to determine whether a substance interferes with the MTT assay, the highest soluble substance concentration is incubated with the MTT solution in the absence of EE. If the MTT assay results in a colour change then the substance must be excluded from further analysis.
- If necessary dilute in isopropanol if the OD is above the limits of the plate reader (e.g. outside the linear range of the plate reader).

Cross reactivity of a substance is determined by the following procedure:

50 µl of the maximum dose of the substance (e.g. 200 mg/ml) is pipetted onto an 8 mm filter paper disk

Treat the filter as if it were an epiCS and perform MTT assay using the filter according to section ENDPOINT MEASUREMENT.

N.B. The test substance should not change the colour of the isopropanol into any other colour. It may not increase OD values more than 2 times the mean OD of background controls. If so, this leads to exclusion of the test substance due to interference with the MTT assay.


Preparation of dose response experiments to determine the EC₅₀ value

In order to determine an EC₅₀ value, substances are tested using two independent epiCS batches in a dose response, starting with the highest soluble concentration (preferably max. 200 mg/ml) and decreasing in 2 fold dilutions until 0.1 mg/ml (or 1st lower concentration hereof) is reached.

Example 1: 200 mg/ml, 150 mg/ml, 100 mg/ml, 50 mg/ml, 25 mg/ml, 12.5 mg/ml, 6.25 mg/ml, 3.125 mg/ml , 1.56 mg/ml, 0.78 mg/ml, 0.39 mg/ml, 0.20 mg/ml, 0.10 mg/ml

Example 2; 150 mg/ml, 75 mg/ml, 37.5 mg/ml, 18.75mg/ml, 9.38 mg/ml, 4.69 mg/ml, 2.34 mg/ml, 1.17 mg/ml, 0.58 mg/ml, 0.29 mg/ml, 0.15 mg/ml, 0.07 mg/ml

Whereas substance dose responses are tested in single fold in each independent experiment, control conditions (unexposed, solvent(s) and positive control), should preferably be tested in duplicate per independent experiment.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 18 of 31	Standard Operating Procedure	

Control conditions and quality criteria for the control conditions

Below a description is given for the different control conditions that are taken along in each experiment. The number of different solvents tested depends on the number of solvents needed for the different substances to be tested within one experiment.

Unexposed epiCS cultures

Unexposed epiCS are cultures that do not receive any treatment during the substance exposure experiment. No filter paper disc, no substance solution and no medium is topically applied to these cultures. In the substance exposure experiments they are taken along as “naked” cultures.

Solvent exposed epiCS cultures

Different solvents, applied to paper filters, have been tested topically on the epiCS: water, PBS, absolute ethanol, 1% DMSO, AOO (4+1) and AOO (1+19)

Quality Criteria for the control conditions

Solvent exposure should not result in more than 30% decrease in relative cell viability compared to unexposed cultures. If the decrease in relative cell viability is > 30% for the solvent exposed cultures, then the epiCS batch does not fulfil the quality criteria required for this assay and the experiment should therefore be excluded from further analysis.

SUBSTANCE EXPOSURE


Unexposed and solvent are only tested in duplicate per batch epiCS, while substances are tested in single.

Preparation of solvents:

- AOO (4+1) - 8 ml acetone + 2 ml olive oil
- AOO (1+19) - 0.5 ml acetone + 9.5 ml olive oil
- 1% DMSO - 50 µl DMSO + 4950 µl epiCS Culture Medium
- PBS – 100% PBS
- Water – sterile dH₂O
- Absolute ethanol

Preparation of testing substance concentrations

- Be aware that the weight / volume of substances is taken into account when dissolving the substances to obtain the Master Stock solution (maximum soluble concentration ≤ 200 mg/ml).


Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 19 of 31	Standard Operating Procedure	

- Substances must be handled with care, as they may be harmful.
- Substances must be freshly prepared before epiCS exposure (to avoid oxidation of the substances, it's recommended to immediately close the tubes firmly and wrap Parafilm around the tubes).
- The use of gloves, lab-coat and mask is recommended.
- For liquid substances, all crystals should be dissolved before pipetting, by leaving the liquid substances at room temperature for up to 1 hour after removal from cooling (if cooling conditions are required for the storage of the substance).

Handling the solvent AOO

It is known that AOO dissolves some plastic lab devices. It is therefore recommended to use polypropylene tubes to make the substance solutions dissolved in AOO. And impregnate the paper filter disks for the substance exposure on glass devices (e.g. glass slides, tissue culture dishes made of glass).

1. Remove old medium from the wells with the inserts and replace it with 1 ml pre-warmed (37°C) epiCS Culture Medium.
2. Place pre-sterilized filter paper disks (8 mm) in a 100 mm tissue culture dish using a sterile pincet (one filter paper disk per epiCS culture). Label the upper side of each filter by marking the filter with a sharp point of a pencil.
3. Impregnate filter paper discs with 25 µl of the testing samples (substance dilutions and solvents (take along duplicates for each testing condition)).
4. Hold the filter paper disk with a pincet and gently tap them minimal 8 times vertically to the plate until the excess of solution slides out from the filter. Be sure to close tissue culture dish immediately to prevent evaporation of solvent in the flow.
5. After impregnating filters, take the epiCS cultures out of the incubator.
6. Quickly (to avoid evaporation of solvent) apply impregnated filters topically to the epiCS stratum corneum according to the plate layout described in Annex 3. Place the filter with the pencil labeled side upwards on the cultures.
7. Incubate exposed cultures for 24 hours (37°C, 5%CO₂, 95% humidity).
8. After 24 hours of incubation, take exposed cultures out of the incubator and gently remove the filter paper disks with a pincet.
9. Perform MTT assay on epiCS to assess cell viability after substance exposure.
10. Transfer the culture medium to 1.5 mL tubes for IL-18 assessment. Store at -20°C or for longer storage at -80°C until ELISA testing.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 20 of 31	Standard Operating Procedure	

ENDPOINT MEASUREMENT:


MTT assay is used to determine the MTT EC₅₀ value (EC₅₀ = effective substance concentration required to reduce cell viability to 50% of the maximum value compared to solvent exposed cultures). IL-18 release is used to characterize the allergenic potential of the tested substance and to give extra potency data (IL-18 SI-2 value).

MTT assay:

MTT [3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide] assay, is based on the ability of a mitochondrial dehydrogenase enzyme from viable cells to cleave the tetrazolium rings of the pale yellow MTT and form dark blue formazan crystals which are largely impermeable to cell membranes. This results in its accumulation within viable cells. Cells will be lysed by the addition of a detergent (isopropanol) resulting in the liberation of the crystals, which are solubilised. The number of surviving cells is directly proportional to the level of the formazan product created. The colour can then be quantified using a scanning spectrophotometer. As MTT is a colorimetric reaction, a scanning spectrophotometer is necessary to quantify the generated colour correlated to cell viability. For incubation steps, an incubator 37 °C, 5% CO₂, 95 % humidity is needed. The MTT assay therefore measures changes in metabolic activity of the cultures, which in turn correlates to changes in viability of the cultures

To perform the MTT assay, follow the steps below:

1. The MTT analysis is performed in a 24-well plate
2. Prepare MTT stock solution in PBS (5 mg/ml): For 50ml stock solution, weigh 250 mg MTT and add 50 ml PBS. Be aware that the PBS that will be used does **NOT** contain Ca²⁺ (calcium) or Mg²⁺ (magnesium). The stock solution can be frozen in aliquots until use. Use MTT working solution at a concentration of 1 mg/ml diluted with PBS.
3. Pipette 300 µl per well of the MTT solution in each well of a 24 well plate.
4. Remove carefully the paper filter discs from epiCS test cultures with a pincet and discard them.
5. Pick up the epiCS test cultures with a pincet and tap the cultures gently on a paper to dry the bottom of the culture. Place the cultures into the 24-well plate previously filled with the MTT solution.
6. Once samples are in contact with MTT solution, incubate for 3 hours in the dark at 37 °C, 5% CO₂, 95% humidity. Make sure that no bubbles are present between the bottom of the cultures and the MTT solution in the well. If so, tap the plate gently to remove the bubbles present underneath the cultures.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 21 of 31	Standard Operating Procedure	

7. It is strongly recommended to take photographs of the 24-wells plate after the incubation of the cultures with the MTT solution, but before the incubation with the isopropanol, in order to document the results also visually, in case any abnormalities may arise from this point on.
8. After the 3 h incubation, pipette 2 ml per well of isopropanol in a new 24-well plate.
9. Pick up the epiCS cultures from the plates with MTT solution. Tap the cultures gently on a paper until the bottom of the culture is dry. Place the cultures into the 24-well plate previously filled with 2 ml isopropanol. Make sure that no bubbles are present between the bottom of the cultures and the isopropanol in the well. If so, tap the plate gently to remove the bubbles present underneath the cultures. Cover the plate with Parafilm to avoid isopropanol evaporation.
10. Incubate the plates with cultures in isopropanol overnight, covered with Parafilm and protected from light or with gently shaking on a plate shaker (~ 100 rpm) for 2 h at room temperature.
11. After the extraction period is complete, pierce the inserts with an injection needle and allow the extract to run into the well from which the insert was taken. Afterwards the insert can be discarded. Before transferring the extract to 96-well plates pipette up and down 3 x until solution is homogenous.
12. Per each tissue transfer 2 × 200 µl aliquots of the blue formazan extracting solution into a 96-well flat bottom microtiter plate. Use isopropanol as blank.
13. Read plates using a spectrophotometer at 540-570 nm wavelength without using a reference filter. Save the plate until step 14 is completed.
14. Calculate cell viability:

In the template files (if provided) in the Excel sheets, formulas and graphs have been incorporated in order to calculate the metabolic activity for each condition.

Absorbance of the wells from the solvent exposed cultures is taken as 100% cell viability. The higher the toxicity of a particular condition, the lower the absorbance value will be. Controls, which are viable have a purple coloured supernatant. The more toxic the condition, the more yellowish the supernatant is and the lower the absorbance value is.


e.g.: Absorbance from control = 0.958

Absorbance from a testing condition = 0.258

0.958 --- 100 % viability

0.258 --- x % viability → cell viability at this specific condition = 26.93%

15. If the absorbance of any of the samples is above 2.0, then samples need to be further diluted and read again with the spectrophotometer. Proceed as follows: for ALL of the samples within a given experiment transfer 100 µl (1:2 dilution) to a new plate and add 100 µl isopropanol. Repeat step 13 and 14 above.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	
Page 22 of 31	Standard Operating Procedure	

ACCEPTANCE CRITERIA

Control conditions and quality criteria for the control conditions

Below a description is given for the different control conditions that are taken along in each experiment. The number of different solvents tested depends on the number of solvents needed for the different substances to be tested within one experiment.

Unexposed epiCS cultures

Unexposed epiCS are cultures that do not receive any treatment during the substance exposure experiment. No filter paper disc, no substance solution and no medium is topically applied to these cultures.

Solvent exposed epiCS cultures

Six different solvents, applied to paper filters, have been tested topically on epiCS: water, PBS, absolute ethanol, 1% DMSO, AOO (4+1) and AOO (1+19).

Quality Criteria for the control conditions

Solvent exposure should not result in more than 30% decrease in relative cell viability compared to unexposed cultures. If the decrease in relative cell viability is > 30% for the solvent exposed cultures, then the epiCS batch does not fulfil the Quality Criteria required for this assay and the experiment should therefore be excluded from further analysis.

DATA ANALYSIS

MTT:

An estimation of the EC50 values can be calculated also as follows:

EC50		Explanation	
conc.	viability	CV<50 (A)	fill in the chemical conc. (in mM, mg/mL or %) with a rel. viability <50%
CV<50 (A)	CV<50 (B)	CV<50 (B)	fill in the rel. cell viability (in %) for the chemical conc. with a rel. cell viability <50%
CV>50 (C)	CV>50 (D)	CV>50 (C)	fill in the chemical conc. (in mM, mg/mL or %) with a rel. viability >50%
EC50 value	50	CV>50 (D)	fill in the rel. cell viability (in %) for the chemical conc. with a rel. cell viability >50%
slope	#WAARDE!	50	this is the 50% relative cell viability
y intercept	#WAARDE!	EC50 value	this is the EC50 value (in mM, mg/mL or %) that induces a rel. cell viability of 50%; the following formula is used to calculate the EC50 value: = (50 - y intercept) / slope
		slope	= (B - D) / (A - C)
		y intercept	= B - (slope * A)

Chemical X	
mg/mL	viability
50,00	28,90
25,00	84,80
40,56	50
slope	-2,24
y intercept	140,70

EC50 value	40,56 mg/mL
------------	-------------


Determining the IL-18 release by ELISA

Expression of IL-18 rationale

IL-18, formerly known as IFN- γ -inducing factor (IGIF), which belongs to the IL-1 cytokine family, is a potent inducer of IFN- γ by activated T cells (Okamura et al., 1995). IL-12 and IL-18 play important roles in the development of T helper type I (Th1) cells and are synergistic in the induction of IFN- γ by T-cells (Okamura et al., 1995). IL-18 has been demonstrated to favour Th-1 type immune response by enhancing the secretion of pro-inflammatory mediators such as TNF- α , IL-8 and IFN- γ , and to play a key proximal role in the induction of allergic contact dermatitis (Cumberbatch et al., 2001). Human keratinocytes constitutively express IL-18 mRNA and protein (Naik et al., 1999). Antonopoulos et al., (2008) have recently demonstrated that IL-18 is a key proximal mediator of LC migration and contact hypersensitivity, acting upstream of IL-1 β and TNF- α , suggestive of a central role in regulation of cutaneous immune response.

IL-18 release is assessed by ELISA. Standards are run in duplicate, samples (0.1 ml) in single.

Results are expressed as pg/ml. The commercially available IL-18 ELISA kit from MBL should be used (#7620, MBL, Japan). Please follow the instructions of the MBL ELISA kit.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 24 of 31	Standard Operating Procedure	

IL-18 Data Analysis

For assessment of allergenicity (YES/ NO): The relative amount of IL-18 protein present in supernatants given in pg/ml. The result is then expressed as fold change, calculated by the following equation:


$$\text{IL-18 SI} = \frac{\text{IL-18 pg/ml in substance-treated EE}}{\text{IL-18 pg/ml in solvent-treated EE}}$$

If the fold increase in intracellular IL-18 is ≥ 5.0 when cell viability $\geq 5\%$ and $\leq 40\%$ (EC_{5-40}) in at least one of the concentrations tested the substance is classified as contact sensitiser (R43).

If the fold increase in intracellular IL-18 is < 5.0 when cell viability $\geq 5\%$ and $\leq 40\%$ (EC_{5-40}) in all concentrations tested the substance is classified as non contact sensitiser.

For a given substance, the same classification must be obtained in both independent experiments. The 5.0 fold increase is indented for at least one of the concentrations tested.

For potency assessment (IL-18 SI-2): the relative amount of IL-18 protein present in supernatants given in pg/ml is determined. The SI-2 is determined using the same spread sheet design as shown above.

Version 1.0 Sep 2014	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 25 of 31	Standard Operating Procedure	

PREDICTION MODEL

A substance is classified as contact allergen if it induces ≥ 5 fold increase in IL-18 release (ELISA) compared to solvent treated epiCS at cell viability $\geq 5\%$ and $\leq 40\%$ (EC_{5-40}) in at least one of the concentrations tested. For sensitiser potency, the EC_{50} value and IL-18 (SI-2) value is determined.


Read out:

For allergenic potential: ≥ 5 fold increase in IL-18 release (ELISA) compared to solvent treated epiCS at cell viability $\geq 5\%$ and $\leq 40\%$ (EC_{5-40}).

For potency assessment:


Primary parameter: Cytotoxicity (MTT assay) expressed as EC_{50} value (EC_{50} = effective substance concentration required to reduce epiCS metabolic activity - corresponding to cell viability - to 50% of the maximum value compared to solvent exposed cultures).

Secondary parameter: IL-18 Substance concentration resulting in : ≥ 2 fold increase in IL-18 release (ELISA) compared to solvent treated epiCS (SI-2).

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	
Page 26 of 31	Standard Operating Procedure	

BIBLIOGRAPHIC REFERENCES/REPORTS

- Antonopoulos, C., Cumberbatch, M., Mee, J.B., Dearman, R.J., Wei, X., Liew, F.Y., Kimber, I., Groves, R.W., 2008. IL-18 is a key proximal mediator of contact hypersensitivity and allergen-induced Langerhans cell migration in murine epidermis. *J. Leukoc. Biol.* 83: 361-367.
- Cavani A, De Pità O, Girolomoni G, 2007. New aspects of the molecular basis of contact allergy. *Curr Opin Allergy Clin Immunol.* 7: 404-408.
- Corsini E, Mitjans M, Galbiati V, Lucchi L, Galli CL, Marinovich M, 2009. Use of IL-18 production in a human keratinocyte cell line to discriminate contact sensitizers from irritants and low molecular weight respiratory allergens. *Toxicol In Vitro.* 23:789-796.
- Corsini E, Galbiati V, Mitjans M, Galli CL, Marinovich M, 2013. NCTC2544 and IL-18 production: A tool for the identification of contact allergens. *Toxicology In Vitro.* 27: 1127-1134.
- Cumberbatch, M., Dearman, R.J., antopoulos, C., Groves, R.W., Kimber, I., 2001. Interleukin-18 induces Langerhans cell migration by a tumor necrosis factor- α and IL-1 β -dependent mechanism. *Immunology* 102: 323-330.
- dos Santos GG, Spiekstra SW, Sampat-Sardjoepersad SC, Reinders J, Scheper RJ, Gibbs S, 2011. A potential in vitro epidermal equivalent assay to determine sensitizer potency. *Toxicol In Vitro.* 25: 347-57.
- Galbiati V, Mitjans M, Lucchi L, Viviani B, Galli CL, Marinovich M, Corsini E, 2011. Further development of the NCTC 2544 IL-18 assay to identify in vitro contact allergens. *Toxicol. In Vitro* 25:724-732.
- Gibbs S, Corsini E, Spiekstra SW Galbiati V, Fuchs HW, DeGeorge G, Troese M, Hayden P, Deng W, Roggen E, 2013. An Epidermal Equivalent Assay For Identification And Ranking Potency Of Contact Sensitizers. *Toxicol. Appl. Pharmacol.* 272: 529-541.
- Naik, S.M., Cannon, G., Burbach, G.J., Singh, S.R., Swerlick, R.A., Wilcox, J.N., Ansel, J.C., Caughman, S.W., 1999. Human keratinocytes constitutively express interleukin-18 and secrete biologically active interleukin-18 after treatment with pro-inflammatory mediators and dinitrochlorobenzene. *J. Invest. Dermatol.* 113: 766-772.
- Okamura, H., Tsutsui, H., Komatsu et al., 1995. Cloning of a new cytokine that induces IFN- γ production by T cells. *Nature* 378: 88-91.
- Spiekstra SW, Dos Santos GG, Scheper RJ, Gibbs S, 2009. Potential method to determine irritant potency in vitro - Comparison of two reconstructed epidermal culture models with different barrier competency." *Toxicol. In Vitro* 23: 349-55.
- Van Och FM, Van Loveren H, Van Wolfswinkel JC, Machielsen AJ, Vandebriel RJ, 2005. Assessment of potency of allergenic activity of low molecular weight compounds based on IL-1 α and IL-18 production by a murine and human keratinocyte cell line. *Toxicology* 210: 95-109.
- OECD TG 429 (2010): OECD Guideline for the testing of substances. Skin Sensitisation: Local Lymph Node Assay.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 27 of 31	Standard Operating Procedure	

ANNEXES

Annex 1

INTERFERENCE TESTS:

Test for interference of substances with MTT as described in the epiCS Skin Irritation INVITTOX Protocol - here adapted to Skin Sensitisation.

Test for interference of substances with MTT endpoint and correction procedures

A test substance may interfere with the MTT endpoint if: a) it is coloured and/or b) able to directly reduce MTT (for possible combination of interactions). The MTT assay is affected only if the test material is present in or on the tissues when the MTT viability test is performed.

Some non-coloured test materials may change into coloured materials in wet or aqueous conditions and thus stain tissues during exposure. Therefore, before exposure, a functional check for this possibility should be performed (Step 1).

Step 1:

Add 25 µl (liquid) of the test substance into 0.3 ml of deionised water. Perform the test in a transparent, preferably glass test tube since plastic test tubes may react with the test articles during the incubation time. Incubate the mixture in the incubator (37±1°C, 5±1 % CO₂, 90±10 % RH) for 60 min. At the end of the exposure time, shake the mixture and evaluate the presence and intensity of the staining (if any). If the solution changes colour significantly, the test substance is presumed to have the potential to stain the tissue. A functional check on viable tissues should be performed (Step 2).

Step 2:

To check the tissue-binding of a coloured test article (or a substance that changes into a coloured substance) the following test should be performed with - we recommend - 3 tissue replicates concurrently tested with each and every test performed with the coloured substance): Expose three tissues with impregnated filter papers. . In parallel, expose a tissue to PBS (negative control). Follow all procedures as described in this protocol in the Method Section except incubate the tissue for 2h in PBS without MTT (37±1°C, 5±1 % CO₂, 90±10 % RH) instead of incubating in PBS containing MTT. After the 2 hr incubation, rinse the tissues and extract the tissues using 2.0 ml of isopropanol and measure the optical density (OD) at 540-570 nm.

Data correction procedure

If the extract from tissues treated by coloured substance (or substance detected in step 1)

has an OD between 5 % and 30 % of the negative control tissue (treated with PBS), the substance should be concurrently tested with each and every test performed with the coloured material using the procedure described above. The OD (unaffected by interference with the coloured test materials) is calculated using following formula:

$$\text{OD} = \text{OD coloured tissue (MTT assay)} - \text{OD coloured tissue (no MTT assay)}$$

Note: If the extract from tissues treated by coloured substance (or substance detected in Step 1) has an OD <5% of the PBS treated control tissue and the tissue viability (determined in MTT assay) is not close to the classification cut-off (50 %), correction of the results is not necessary.

If the OD of extract from the tissue treated by coloured substance (or substance detected in Step 1) is > 30 % of the PBS treated control tissue, additional steps and expert judgment must be performed to determine if the test substance must be considered as incompatible with the test.

Step 3:

All test materials (including those already evaluated in Step 1 and Step 2) should be further evaluated for their potential to interfere with MTT assay. To test if a material directly reduces MTT, add 25 µl test substance to 1 ml PBS containing MTT and incubate in the incubator (37±1°C, 5±1 % CO₂, 90±10 % RH) for 60±2 min. As negative control use 1 ml PBS containing MTT. If the MTT solution turns blue/purple, the test substance reduces MTT and additional functional check (Step 4) must be performed.





Figure 1: Example of test for direct MTT reduction ability (Step 3). Test substances (2) (3) and (6) have directly reduced MTT. In these cases, Step 4 must be performed.

Version 1.2 Nov. 2015	epiCS[®] In Vitro Skin Sensitisation INVITTOX Protocol	 CELLSYSTEMS[®]
Page 29 of 31	Standard Operating Procedure	

Step 4:

The procedure employs freeze-killed tissues that possess no metabolic activity but absorb and bind the test substance similar to viable tissues.

Prepare freeze-killed tissues:

- a) transfer epiCS in 24-well plates without cell culture medium
- b) freeze tissue at -20°C (or -80°C) for at least 12 h (3 tissues / MTT-interacting test substance)
- c) thaw tissues and repeat step b) twice
- d) thaw tissue 1 h (± 10 min) before use
- e) keep at RT

Correction of non specific MTT reduction should be performed with 3 tissue replicates in a single run, independently of the number of independent tests performed with the direct MTT-reducing substance:

Each MTT reducing substance is applied to three freeze-killed tissue replicates. In addition, three freeze killed tissues are left untreated (Note: The untreated killed controls will show a small amount of MTT reduction due to residual reducing enzymes within the killed tissue). The entire assay protocol is performed with the frozen tissues in parallel to the viable epiCS tissues. Data are then corrected as follows:

Data correction procedure

True viability = Viability of treated tissue – MTT conversion by substance

= OD tvt – OD kt

OD kt = (mean OD tkt – mean OD ukt)

tvt = treated viable tissue

kt = killed tissues

tkt = treated killed tissue

ukt = untreated killed tissue (NC treated tissue)

If the interference by the test substance is greater than 30 % of the negative control value, additional steps must be taken into account or the test substance may be considered incompatible with this test system

If the interference by the test substance is ≤ 30 % of the negative control value, the net OD of the test substance treated killed control may be subtracted from the mean OD of the test substance treated viable tissues to obtain the true amount of MTT reduction that reflects metabolic conversion only.

Annex 2

Dissolving substances in the solvent acetone:olive oil 4+1 (AOO)

Dose Response for AOO soluble substances (100% final AOO)

→ In the dose response 12x 2-fold serial dilutions (1x, 2x, 4x, 8x, 16x etc) are tested for each substance.

Steps to be taken


1. Make AOO (4+1): 8 ml acetone + 2 ml olive oil
2. Dilutions are made in AOO as indicated in the table below starting with the Master Stock concentration
3. Impregnate the filter paper discs with the substance dilutions and expose and analyse the cultures according to the protocol described.

Based upon the results from the dose response experiments an EC₅₀ value for each substance might be determined.

Dilution to be tested	µl of previous dilution	µl of AOO (4+1)
1x dilution	Master Stock concentration in 100% AOO (e.g. 200 mg/ml)	--
2x dilution	300 µl (1x dilution)	300 µl
4x dilution	300 µl (2x dilution)	300 µl
8x dilution	300 µl (4x dilution)	300 µl
16x dilution	300 µl (8x dilution)	300 µl
32x dilution	300 µl (16x dilution)	300 µl
64x dilution	300 µl (32x dilution)	300 µl
128x dilution	300 µl (64x dilution)	300 µl
256x dilution	300 µl (128x dilution)	300 µl
512x dilution	300 µl (256x dilution)	300 µl
1024x dilution	300 µl (512x dilution)	300 µl
2048x dilution	300 µl (1024x dilution)	300 µl

Annex 3

Control plate



Test chemical dilution plates

